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
WATER DEMINERALIZER

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USER MANUAL

REVISIONS

The following table lists subsequent editions/revisions of the manual.
The "Description" field brief explains the subject of the latest revision.

Edit.	Rev.	Date	Description
1	0	09-2005	First issue
1	1	10-2006	Modified resin filter assembly and maintenance procedures

TABLE OF CONTENTS

INTRODUCTION	1
APPLICABLE EUROPEAN DIRECTIVES	1
INTENDED USE	1
PURPOSE OF THE MANUAL	2
GENERAL DIRECTIONS	2
PACKAGE CONTENT	3
SIZE AND WEIGHT	3
HANDLING THE PRODUCT	3
CONTENT DESCRIPTION	3
PRODUCT INFORMATION.....	4
INTRODUCTION	4
GENERAL FEATURES	4
FRONT PANEL	5
REAR PANEL	5
FUNCTIONING	5
INSTALLATION	6
INTRODUCTION	6
POWER SUPPLY	6
SIZE OF EMBEDDING CHAMBER	6
GENERAL INSTALLATION PRECAUTIONS	7
ELECTRIC CONNECTIONS	7
INSTALLATION ONTO A WALL	8
HYDRAULIC CONNECTIONS	9
Initial installation	9
Purging the internal hydraulic circuit (osmotic membrane)	10
Case 1: Direct connection to a sterilizer series Millenium	11
Case 2: Connection to an external feeding tank (optional)	12
Case 3: Connection to an electric-type drawing tap (option)	13

FIRST STARTUP	15
MILLDROP CONNECTED TO THE STERILIZER SERIES MILLENNIUM.....	15
MILLDROP CONNECTED TO THE FEEDING TANK	15
MILLDROP CONNECTED TO THE MANUAL DRAWING TAP	15
STAND-BY STATUS	16
APPENDIX A – TECHNICAL CHARACTERISTICS	17
SUMMARY TABLE	17
CHARACTERISTICS OF THE DEMINERALISED WATER.....	18
APPENDIX B - MAINTENANCE	19
INTRODUCTION	19
ORDINARY MAINTENANCE PROGRAMME	19
GENERAL MAINTENANCE PROCEDURES	19
REPLACING THE FILTERS	20
Resin filter cartridge	20
Osmotic membrane.....	21
Active carbon filter	22
REPLACING THE FLOW RESTRICTOR	23
APPENDIX C – TROUBLESHOOTING.....	24
INTRODUCTION	24
ANALYSIS AND RESOLUTION OF PROBLEMS	24
APPENDIX D - DIAGRAMS	26
ELECTRIC DIAGRAM	26
Electronic board type 1	26
Electronic board type 2	27
HYDRAULIC DIAGRAM	28
APPENDIX G – DECLARATION OF CONFORMITY	29
APPENDIX H - NOTES	30
APPENDIX Z – TECHNICAL ASSISTANCE.....	31

INTRODUCTION

Dear Client

Thank you for choosing this product from M.O.COM. Srl. We hope that you will find it completely satisfactory.

This manual describes all procedures for the correct installation and use of the equipment.

Should you have any questions or suggestions concerning this product or our support, do not hesitate to contact us.

Symbols

NOTE

PAY SPECIAL ATTENTION; NOTES FOR PROPER INSTALLATION AND MAINTENANCE PROCEDURES AND USE METHODS.

WARNING

THIS SYMBOL INDICATES A POTENTIAL DANGER OF INJURY. FOLLOW THE PROCEDURES DESCRIBED IN THE MANUAL TO AVOID INJURING THE USER AND/OR OTHERS.

DANGER

THIS SYMBOL SIGNALS A POSSIBLE DANGER FOR PEOPLE. FOLLOW THE SUGGESTED PROCEDURE IN ORDER TO PREVENT POSSIBLE DAMAGE TO USERS AND/OR OTHERS.

THE MATERIAL THE DEMINERALIZER IS COMPOSED OF MUST BE DISPOSED ACCORDING TO THE DIRECTIVE 2002/96/CEE

ENGLISH

APPLICABLE EUROPEAN DIRECTIVES

The product here described is made according to the highest safety standards and does not present any danger for the operator, provided it is used according to the directions reported hereafter.

The product complies with the following applicable European Directives:

73/23/EEC (and subsequent amendments), for reconciling Member State legislations regarding low tension apparatuses;

89/336/EEC (and subsequent amendments), for reconciling Member State legislations regarding electromagnetic compatibility.

INTENDED USE

The equipment is meant for professional use and can be employed only the purpose it was designed for, that is the demineralization of the public water, and its use for the water steam sterilizer feeding.

Important notes

WARNING

THE DEVICE MUST ONLY BE USED BY QUALIFIED PERSONNEL. IT MAY NOT BE USED OR HANDLED BY INEXPERT AND/OR UNAUTHORIZED PERSONNEL FOR ANY REASON.

NOTES

THE INFORMATION INCLUDED IN THIS MANUAL IS SUBJECT TO CHANGE WITHOUT NOTICE.

M.O.COM. SRL SHALL NOT BE RESPONSIBLE FOR DIRECT, INDIRECT, ACCIDENTAL OR CONSEQUENTIAL DAMAGE OR OTHER DAMAGES RELATED TO THE SUPPLY/USE OF THIS INFORMATION.

THIS DOCUMENT MAY NOT BE REPRODUCED, ADAPTED OR TRANSLATED, IN WHOLE OR IN PART, WITHOUT THE PRIOR, WRITTEN AUTHORIZATION OF **M.O.COM. SRL**

mocom IS REGISTERED TRADEMARK OF **M.O.COM. SRL**.

PURPOSE OF THE MANUAL

The purpose of this manual is to provide directions for:

- becoming familiar with the product;
- its appropriate installation;
- its secure and efficient employment.

The appendixes provide information regarding:

- technical characteristics;
- maintenance procedures;
- analysis and solution of problems;
- further miscellaneous documentation.

GENERAL DIRECTIONS

The product must **always** be used in compliance with the procedures set out in this manual and never for purposes that are different from the foreseen ones.

WARNING

THE USER IS RESPONSIBLE FOR ALL LEGAL OBLIGATIONS LINKED TO THE PRODUCT'S INSTALLATION AND USE. IF THE PRODUCT IS NOT INSTALLED OR USED PROPERLY, OR CORRECT MAINTENANCE/ SERVICING IS NOT CARRIED OUT, IT SHALL NOT BE POSSIBLE TO DEEM THE MANUFACTURER RESPONSIBLE FOR ANY POSSIBLE BREAKAGE, MALFUNCTIONING, DAMAGES TO THINGS, OR INJURIES TO PERSONS.

Please observe the following precautions in order to avoid injury or property damage; in particular:



- Before connecting the power supply cable, **ensure** that nominal mains voltage is equal to the one indicated on the label. Voltages differing from indicated ones, may damage the apparatus's electric parts.
- Be sure the electric plant is provided with a **grounding** conforming to laws and/or rules in force.
- When getting started, carry out the "hydraulic circuit purging" so that possible processing residuals will be expelled (see "Initial installation").

DANGER

ALWAYS REMOVE THE MAINS CABLE FROM THE POWER OUTLET, BEFORE CARRYING OUT ANY KIND OF MAINTENANCE. SHOULD THAT NOT BE POSSIBLE, OR WHENEVER THE EXTERNAL BREAKER IS FAR AWAY OR, AT ANY RATE, NOT VISIBLE TO THE MAINTAINER, PLACE A WORK IN PROGRESS SIGN ON THE EXTERNAL BREAKER AFTER TURNING IT OFF .

- **Do not** remove any label or nameplate from the device; should you need any new ones, apply for them.
- **Do not** use corrosive products, acids, steel wool or brushes for cleaning the equipment.
- **Do not** wash the equipment with direct water jets.
- Use **only original spare parts**.



PACKAGE CONTENT

SIZE AND WEIGHT


Package:

Overall weight about 18 Kg

Once the package is opened, check that:

- the supply matches the specifications of the order (see the accompanying document);
- there is no obvious product damage.

ENGLISH

HANDLING THE PRODUCT

The product must be handled according to directions provided on the packaging.

We recommend to transport and store the device at temperatures not lower than 5°C. Prolonged exposure to low temperature can damage the product.


Keep the original packaging and use it whenever the device is to be transported. The use of different packaging could damage the product during shipment.

Keep the plugs mounted on the demineralizer's pipe fittings (to be reapplied in case of subsequent shipment).

CONTENT DESCRIPTION

The package contains the following items:

- Demineralizer equipment (Ref. 1);
- Kit for the installation of the draining pipe (Ref. 2);
- Draining pipe (red) - Ref. 3;
- Feeding pipe (blue), complete with gaskets – Ref. 4;
- Reinforced silicone pipe (Ref. 5) for connecting MILLDROP to sterilizer, external tank or tap;
- Power supply cable – Ref. 6;
- Kit (Ref. 7) containing: start/stop signal cable, pipe union, nr. 4 clamps;
- Operator's Manual and Warranty Certificate – Ref. 8.

PRODUCT INFORMATION

INTRODUCTION

"MILLDROP" is **M.O.COM.**'s proposal in the field of water demineralization for water steam sterilizer.

It is a sophisticated apparatus but easy to use, and adaptable to the user requirements by means of different installation possibilities.

Thanks to the reverse osmosis membrane and resin filter, the equipment "transforms" the water coming from water works into qualitatively excellent demineralized water.

"MILLDROP" – which is considerably easy to use and small in size - represents the ideal equipment to be put next to any water steam sterilizer.

GENERAL FEATURES

"MILLDROP" is a water demineralizer based on the reverse osmosis principle.

The equipment includes:

- active carbon filter, having a duration of about 800 liters of demineralized water produced;
- two resin filters having a duration of about 500 liters of demineralized water produced,
- osmotic membrane which, if regularly checked, can last more than 800 liters of demineralized water produced.

NOTE

THE LIFETIME OF FILTERS AND OSMOTIC MEMBRANE DEPENDS ON THE QUALITY OF THE WATER FEEDING THE MILLDROP.

The limit values for the feeding water are reported in the following table:

Element/Feature	Value
Turbidity	< 1 NTU
Choline	< 0.1 ppm
Ph	2-11
Salt density	< 5 SDI
Specific conductivity at 20°C	< 1000 µS/cm
Iron	< 0.1 ppm
Temperature	< 35°C
Total hardness	< 30°F

NOTE

FOR A SPECIFIC HARDNESS HIGHER THAN 25°F, THE USE OF A RESIN SWEETENER AND SALTS REGENERATION IS RECOMMENDED ; THIS IS REQUIRED IF THE HARDNESS IS HIGHER THAN 30°F.

The equipment operation is controlled by an electronic board which also monitors any alarm (leak and lack of water, quality of the output water, etc.).


"MILLDROP" may be installed on a work plane, hung to a wall or set inside the washbasin piece of furniture.

"MILLDROP" offers different installation possibilities:


- direct connection to the sterilizer series Millennium (provided with automatic filling feature and start/stop interface);
- direct connection to an external tank with automatic water filling feature (possibility to connect up to 3 sterilizers);
- direct connection to a drawing tap (provided with electric button), to fill manually little tanks.

For extra details please refer to chapter "Hydraulic connections".

FRONT PANEL

REAR PANEL

FUNCTIONING

The MILLDROP operation is very easy and automatic in any use.

- 1) Directly connected to a Millennium series sterilizer, it operates at the start of the sterilizer's water feeding pump, and until the pump stops.
- 2) Connected to the external tank, it starts working as the water level in the tank goes down the maximum threshold, and stops as it is again reached.
If continually working for a time exceeding 20 minutes, MILLDROP shall stop automatically (see Note).
- 3) Connected to the drawing tap (provided with electric device), it starts working as you press the tap's button and stops pressing it again.
If continually working for a time exceeding 20 minutes, MILLDROP shall stop automatically (see Note).

NOTE

TO RESTART THE DEMINERALIZER FROM A STOP STATUS, TURN IT OFF AND ON AGAIN.

INSTALLATION

INTRODUCTION

The correct and careful installation is the first and fundamental step in achieving the good demineralizer operation, long life and complete utilization of its features. Moreover, this precaution will avoid the danger for physical injury or property damage, not to mention equipment malfunctions and faults. So, please follow the instructions in this chapter scrupulously.

NOTE

THE TECHNICAL SERVICE (SEE APPENDIX Z) IS AT YOUR DISPOSAL FOR ANY DOUBT YOU MAY HAVE OR FURTHER INFORMATION YOU MAY REQUIRE.

THE DEMINERALIZER HAS PASSED THE REQUIRED TESTS BEFORE BEING DELIVERED. SO IT DOES NOT REQUIRE ANY ADDITIONAL CALIBRATION AT THE GETTING IN SERVICE.

POWER SUPPLY

The electrical system the demineralizer will be connected to, must be suitably dimensioned basing on the power requirements of the device. See the nameplate located on the backside of the equipment.

WARNING


VERIFY THE ELECTRIC PLANT'S CONFORMITY, AND THAT IT IS GROUNDING PROVIDED.

SIZE OF EMBEDDING CHAMBER

When installing the demineralizer inside a cabinet, adequate space all around the device must be foreseen in order to provide an effective ventilation, as well as an opening in the back to allow the passage of the power cord and piping.

The built-in compartment should have the following minimum dimensions:

- Height 280 mm
- Width 270 mm
(free space on the ventilation grid side)
- Depth 480 mm
(plus a space for rear connections and front panel view)

WARNING

COMPARTMENT DIMENSIONS LESS THAN THOSE SHOWN MAY COMPROMISE THE CORRECT AIR CIRCULATION AROUND THE DEVICE AND AN ADEQUATE COOLING, CAUSING CONSEQUENT PERFORMANCE DETERIORATION AND/OR POSSIBLE DAMAGE.

NOTE

IF THE MAINS SWITCH REMAINS INACCESSIBLE, USE A MAINS SOCKET INCORPORATING AN ON/OFF SWITCH.

REFER TO APPENDIX A (TECHNICAL CHARACTERISTICS) FOR COMPLETE TECHNICAL DATA.

GENERAL INSTALLATION PRECAUTIONS

Obey the following warnings for the correct operation of the device and/or avoid **dangerous situations**:

- Install the demineralizer on a **flat surface**;
- Leave adequate space for ventilation (at least 5 cm on the of the ventilation grid side).
- Do not install near cisterns, wash-basins, etc. to avoid direct contacts with water or liquids, possible short circuits and/or potentially dangerous situations for the operator;
- Do not install in excessively humid or poorly ventilated places;
- Do not install in places with **gas** or inflammable and/or explosive **vapors**;
- Install the unit so that the power cord is not bent or crushed. It must run freely all the way to the mains socket.
- Install the unit so that any external fill/drain piping is not bent or crushed. The length of the reinforced silicone tube connecting MILLDROP to the sterilizer or tank, **shall not exceed 2 meters**, otherwise performance quality could fall or, what is worse, the demineralizer itself could not operate properly.

ENGLISH

ELECTRIC CONNECTIONS

The demineralizer must be connected to a mains socket of adequate power for the absorption required, and provided with ground connection in compliance with current laws and/or standards.

The mains socket must be suitably protected by a circuit breaker having the following characteristics:

- Nominal current **5 A**
- Differential current **0,03 A**

NOTE

THE MANUFACTURER WILL NOT BE LIABLE FOR DAMAGES CAUSED BY INSTALLING THE UNIT ON INADEQUATE ELECTRICAL SYSTEM AND/OR NOT GROUNDED.

If you need to replace the power cord plug, use a type of equal characteristics or, at any rate, adequate to the equipment power requirements. The user is entirely responsible for the choice and replacement of the plug.

DANGER

CONNECT THE POWER CORD DIRECTLY TO THE SOCKET.
DO NOT USE EXTENSIONS CORDS, ADAPTERS OR OTHER ACCESSORIES.

INSTALLATION ONTO A WALL

MILLDROP can be placed on a working plane, mounted on the Millrack system (see the relating User Manual) or installed on a wall (see Fig. A-B-C). In this last case, proceed as follows:

1. Remove the screws (1) and the lid (A). Figure A
2. Remove the screws (2-3) and the side panel (B). Figure A
3. Position the panel (B) on the wall where MILLDROP has to be installed and mark the wall by a pencil in the mounting slots (C). Figure B

NOTE

THE DEMINERALIZER MUST BE INSTALLED PERFECTLY HORIZONTAL, OTHERWISE THE ANTI-FLOODING SYSTEM COULD NOT PROPERLY FUNCTION. ON CHOOSING THE INSTALLATION POINT, CONSIDER THE ROOM NECESSARY FOR REAR PIPING, FOR THE FRONT PANEL FREE VISION AND FOR THE ADEQUATE VENTILATION (SEE PARAGRAPH GENERAL INSTALLATION PRECAUTIONS).

4. Bore the holes corresponding the marks just made. Figure C.

NOTE

BEFORE DRILLING, BE SURE THAT NO PIPING OR ELECTRIC CABLES ARE INSIDE THE WALL.

5. Reassemble the side panel (B), fasten the screws (2 and 3).
6. Insert the screw anchors (4) in the holes, the screws (5) in the panel's slots and screw down. (Note: items 4 and 5 are not supplied) – see Figure C.
7. Reassemble the lid (A) and fasten the screws (1).


Figure A


Figure B


Figure C

HYDRAULIC CONNECTIONS

Initial installation

See Figures D, E:

1. Use the kit supplied and install the draining bracket onto the washbasin's siphon where you wish to install the MILLDROP draining pipe; carry out an hole Ø 6 mm in the siphon across the bracket's fitting.
2. Insert the red pipe in the ring nut , plug it in the bracket's fitting, and fasten the ring nut.


Figure D

3. Cut the red pipe long enough to reach the draining connection B on the rear panel of the demineralizer, and plug-in the pipe.

NOTE

THE PIPING MUST NOT BE BENT, SQUASHED OR HAVE A PARTICULARLY TORTUOUS ROUTE, OTHERWISE PERFORMANCE WILL BE QUALITATIVELY VERY LOW OR, WHAT IS WORSE, THE DEMINERALIZER ITSELF SHOULD MALFUNCTION.

4. Install a tap (not supplied) provided with male connector 3/4" on the water outlet and near the demineralizer. Check the pressure value of the water works (see Appendix A).
5. Connect the feeding pipe (blue) between the tap outlet and the water inlet A on the rear panel of the demineralizer (check the presence of gaskets), fasten both the pipe's ring nuts.
6. Open the tap and check for no leaks.
7. Remove the plug from the outlet C on the rear panel of the demineralizer and plug-in the pipe fitting.
8. Mount the reinforced silicone pipe on the pipe fitting, tight by clamps.
9. Place the opposite end of the silicone pipe inside a wash-basin.


Figure E

Purging the internal hydraulic circuit (osmotic membrane)

See Figure F and G:

1. Remove the screws and the lid (see Figure A) to access the electronic board (see Figure F)


Electronic board type 1:

- Remove the protection panel of the electronic board, unplug the signal terminal board.
- Remove the blue wires (coming from the relay board for the sterilizer automatic filling interface) from pins 5-6 of the terminal board.
- Make a "jumper" between pins 5-6 and fit the terminal board on the electronic board.

Electronic board type 2:

- Remove (if present) the wires from pins 3-4 of the terminal board, and make a "jumper".

Electronic board type 1

Electronic board type 2


Figure F

2. Disconnect the pipe from the input resin filter fitting, extend it by using the silicone pipe supplied and take its end in a washbasin (figure G).


Figure G

**Case 1:
Direct connection to
a sterilizer series
Millenium**

3. Connect the mains cable and switch on the demineralizer.
4. Now the water will start to flow from the silicone tube and from the red pipe connected to fitting B on the demineralizer back panel; let to drain for about 10 minutes in order to expel any possible processing residue.
5. Turn off the demineralizer and disconnect the mains cable.
6. Restore the internal hydraulic and electric connections and replace the cover.

Continue the installation by carrying out the specific hydraulic connections according to the use required.

See Figure H:

1. Cut the silicone pipe (connected to the fitting C on the demineralizer rear panel) long enough to reach the pipe fitting installed on the automatic feeding point of the sterilizer, and tight by clamps.

NOTE

THE PIPING MUST NOT BE BENT, SQUASHED OR HAVE A PARTICULARLY TORTUOUS ROUTE. THE LENGTH OF THE SILICONE TUBE MUST NOT EXCEED 2 METRES, OTHERWISE PERFORMANCE WILL BE QUALITATIVELY VERY LOW OR, WHAT IS WORSE, THE DEMINERALIZER ITSELF SHOULD MALFUNCTION.

2. Connect the supplied start/stop signal cable between the connector S on the demineralizer rear panel and the female jack on the sterilizer.


Figure H

**Case 2:
Connection to an
external feeding tank
(optional)**

See Figure I and L:

1. Access the electronic board:


Electronic board type 1:

- Unplug the signal terminal board and disconnect from pins 5-6 the blue wires coming from the relay board; insulate the wires by tape.
- Disconnect the red wires from the relay board's terminals and connect them to pins 5-6 of the signal terminal board; tighten the screws and fit the terminal board on the electronic board.

Electronic board type 2:

- Disconnect the read wires from the relay board, connect them to pins 3-4 of the terminals board and tighten the screws.

Electronic board type 1

Electronic board type 2


Figure I

2. Plug the supplied pipe fitting in the inlet (label "MILLDROP") of the tank.
3. Cut the output silicone pipe (connected to the outlet C on the demineralizer rear panel) long enough to reach the pipe fitting installed on the tank's inlet (label "MILLDROP"), and tight by clamps.

NOTE

THE PIPING MUST NOT BE BENT, SQUASHED OR HAVE A PARTICULARLY TORTUOUS ROUTE. THE LENGTH OF THE SILICONE TUBE MUST NOT EXCEED 2 METRES, OTHERWISE PERFORMANCE WILL BE QUALITATIVELY VERY LOW OR, WHAT IS WORSE, THE DEMINERALIZER ITSELF SHOULD MALFUNCTION.

4. Plug a supplied pipe fitting in one of the three outlets available on the tank (label M1, M2 or M3), fit a silicone pipe and tight by clamps.
5. Cut the silicone pipe long enough to reach the pipe fitting installed on the sterilizer automatic feeding point, and tight by clamps.
6. Connect the start/stop signal cable (provided with the demineralizer) between the connector S on the demineralizer rear panel and the female jack of the tank's level sensor cabling.
7. The tank is provided with a tap for drawing manually the permeation product. If the user wishes to use also this feature, connect the supplied pipe fitting to this drawing tap.


Figure L

**Case 3:
Connection to an
electric-type drawing
tap (option)**

See Figure M and N:

1. Access the electronic board and modify the wire connection (see "Connection to an external feeding tank").
2. Change also the jumper setting as follows:


Figure M

3. Reassemble the board protection panel and the lid of the demineralizer.

4. Install the drawing tap (with electric button) where desired (anyhow, near the demineralizer).
5. Plug the supplied pipe fitting on the drawing tap's outlet, cut the silicone pipe installed on the outlet C of the demineralizer so that it is long enough to reach the pipe fitting of the tap, and tight by clamps.

NOTE

THE PIPING MUST NOT BE BENT, SQUASHED OR HAVE A PARTICULARLY TORTUOUS ROUTE. THE SILICONE TUBE'S LENGTH MUST NOT EXCEED 2 METRES, OTHERWISE PERFORMANCE WILL BE QUALITATIVELY VERY LOW OR, WHAT IS WORSE, THE DEMINERALIZER ITSELF SHOULD MALFUNCTION.

6. Connect the start/stop signal cable (provided with the demineralizer) between the connector S on the demineralizer rear panel and the button connector of the tap.


Figure N

FIRST STARTUP

MILLDROP CONNECTED TO THE STERILIZER SERIES MILLENNIUM

Once the MILLDROP has been properly installed, proceed with the following additional startup operations and checks according with the type of use.

Turn on the demineralizer and then the sterilizer.

Ended the initial test of the sterilizer, enter the SETUP menu, select "Filling Options" and the mode "AUTOMATIC" (see the Operating Manual of the sterilizer).

Check that demineralizer starts to work as the sterilizer's automatic feeding pump starts, and stops at the same time as the pump.

NOTE

THE EQUIPMENT WILL STOP AUTOMATICALLY IF THE DELIVERY OF WATER LASTS UNINTERRUPTEDLY FOR MORE THAN 20 MINUTES. TO RESTART THE MILLDROP OPERATION TURN OFF AND ON THE UNIT AGAIN.

MILLDROP CONNECTED TO THE FEEDING TANK

The demineralizer starts immediately as turned on.

Reached the maximum level in the tank (signaled through the level sensor), the demineralizer stops.

Turn on the sterilizer.

Ended the sterilizer initial test, enter the SETUP menu, select "Filling Options" and the mode "AUTOMATIC" (see the Operating Manual of the sterilizer).

If the water level into the sterilizer's internal tank is lower the minimum, the automatic feeding pump starts.

In this case, the demineralizer will start a few seconds after the start of the sterilizer's pump (i.e. as the tank's float, by lowering, supplies the start/stop command to the demineralizer), and will stop as the tank's float detects again the maximum threshold.

NOTE

THE EQUIPMENT WILL STOP AUTOMATICALLY IF THE DELIVERY OF WATER LASTS UNINTERRUPTEDLY FOR MORE THAN 20 MINUTES. TO RESTART THE MILLDROP OPERATION TURN OFF AND ON THE UNIT AGAIN.

In case of manual water drawing also from the tank's tap, ensure that the water remaining in the tank is sufficient to warrant the need of sterilizer(s) connected to the tank, considering also the time necessary to restore a sufficient level in the tank. In substance:

- one sterilizer connected, the remaining water must be at least **2 liters** .
- two sterilizers connected, the remaining water must be at least **4.5 liters**.
- three sterilizers connected, the remaining water must be at least **7 liters**.

MILLDROP CONNECTED TO THE MANUAL DRAWING TAP

Press the tap's button to start the water delivery.

Press again the tap's button to stop the water delivery.

NOTE

THE EQUIPMENT WILL STOP AUTOMATICALLY IF THE DELIVERY OF WATER LASTS UNINTERRUPTEDLY FOR MORE THAN 20 MINUTES. TO RESTART THE MILLDROP OPERATION TURN OFF AND ON THE UNIT AGAIN.

STAND-BY STATUS

Over the first startup operations, the demineralizer sets in stand-by status.

In such a condition, the turning on of any LED will have no influence on the demineralizer subsequent starting, except the alarm "Internal leak" that, if active, stops the demineralizer and prevents it from operating.

In this case, turn off the demineralizer and call the Technical Service.

NOTE

THE TECHNICAL SERVICE (SEE APPENDIX Z) IS AT YOUR DISPOSAL FOR ANY DOUBT YOU MAY HAVE OR FURTHER INFORMATION YOU MAY REQUIRE.

SUMMARY TABLE

Device	REVERSE OSMOSIS WATER DEMINERALIZER
Model	MILLDROP
Manufacturer	C.M.F. S.r.l. Via E.Ellero, 9 33080 Cusano di Zoppola (PN [=Province of Pordenone]) - ITALY for M.O.COM. S.r.l. Via delle Azalee, 1 20090 Buccinasco (MI [=Province of Milan]) – ITALIA
Power supply voltage	220 – 240 V ~
Mains frequency	50 Hz
Mains fuses	2x F 3,15A 250V
Fuse on the electronic board	T 3,15A 250V
External size (LxHxD)	220 x 235 x 415 mm (rear connections excluded)
Nominal power	130 W, 2A max
Insulation class	Class I
Utilization	Internal use
Environment operating conditions	Temperature: +15°C ÷ +35°C Relative humidity: 90% non-condensing
Net weight	14 Kg
Inlet flow rate	minimum 100 lt/h
Outlet flow rate	25 lt/h
Inlet pressure	1 ÷ 4 bar
Working pressure	Max 10 bar (at pump outlet)
Tap water temperature	+10°C ÷ + 35°C
Resin filter cartridge life (average value)	about 500 l (depending on the quality of the inlet water)
Osmotic membrane and active carbon filter life (average value)	about 800 l (depending on the quality of the inlet water)

CHARACTERISTICS OF THE DEMINERALISED WATER**FILTERING CHARACTERISTICS**

The following table shows the demineralizer's average filtering characteristics, compared with values stated by the standards EN 285 and EN 13060.

Element/Characteristic	Unit of measure	Value at the inlet	Value at the outlet	Value according to EN 13060/EN 285
Ph value	-	7,6	6,3	5 - 7
Conductivity at 20°C	µS/cm	435	1,8	<15
Dry residue	mg/l	345	2	<10
Cadmium	mg/l Cd	<0.01	<0.003	<0.005
Iron	mg/l Fe	0.32	0.03	<0.2
Manganese	mg/l Mn	0.03	<0.003	<0.1
Nickel	mg/l Ni	<0.01	<0.003	<0.1
Lead	mg/l Pb	<0.05	<0.02	<0.05
Copper	mg/l Cu	0.02	0.006	<0.1
Zinc	mg/l Zn	0.56	0.02	<0.1
Chlorides	mg/l Cr	3.7	<0.05	<2
Phosphates	mg/l P	0.01	<0.01	<0.5
Silicon oxide	mg/l SiO ₂	5,7	<0.10	<1
Appearance	-	colorless, transparent, without sediments		
Hardness	mmol/l	2.57	<0.005	<0.02


NOTE

THE VALUES SHOWN IN THE COLUMN "OUTLET WATER" ARE INDICATIVE AND REFERRED TO VALUES OF THE COLUMN "INLET WATER".

INTRODUCTION

In addition to the procedures for the correct use, it needs the user performs an ordinary maintenance in order to guarantee safe and efficient operation over the entire life of the equipment.

The ordinary maintenance hereafter described consists in easy manual activities and preventive maintenance interventions.

**ORDINARY
MAINTENANCE
PROGRAMME**

The table summarizes the maintenance required to maintain the demineralizer operating at its peak efficiency.

WEEKLY	Clean and disinfect the external surfaces
YEARLY	Replace the active carbon filter (or earlier if the limit of 800 litres has been reached)
	Replace the resin filter (or earlier if the limit of 500 litres has been reached)
	Replace the osmotic membrane (or earlier if the limit of 800 litres has been reached)

DANGER

ALWAYS REMOVE THE MAINS CABLE FROM THE POWER OUTLET, BEFORE CARRYING OUT ANY KIND OF MAINTENANCE.

SHOULD THAT NOT BE POSSIBLE, OR WHENEVER THE EXTERNAL BREAKER IS FAR AWAY OR, AT ANY RATE, NOT VISIBLE TO THE MAINTAINER, PLACE A WORK IN PROGRESS SIGN ON THE EXTERNAL BREAKER AFTER TURNING IT OFF .

**GENERAL
MAINTENANCE
PROCEDURES**

With reference to the above table, let's take a summary look at the various maintenance to be performed:

- Clean all the external surfaces using a clean cotton cloth dampened with water, and , possibly, the addition of a neutral detergent. Dry the surfaces and remove any residue before using the equipment.
- To disinfect the external surfaces, use denatured alcohol, as well as detergents containing a minimum percentage of sodium hypochlorite (or equivalent substances).
- Yearly or whenever necessary, replace the exhausted filters with new ones having the same characteristics. For this maintenance, call the Technical Service (see Appendix Z).

NOTE

SHOULD THE DEMINERALIZER BE FORESEEN UNUSED FOR A TIME LONGER THAN 20 DAYS, WE SUGGEST TO EMPTY THE EXTERNAL FEEDING TANK, IF CONNECTED.

REPLACING THE FILTERS

Resin filter cartridge

For filters replacement follow the procedures reported below.

See Figures O, P:

1. Turn off MILLDROP, remove the mains cable and disconnect pipes and cable from the use side (sterilizer, feeding tank or drawing tap).
2. Turn off the tap feeding the water to MILLDROP.
3. Remove the screws (1) and take off the lid (A).
4. Remove the screws (2) and take off the side panel (B) (if possible).
5. Lift and extract the filters assembly (C) from its bearings, disconnect pipes Ø8 (3 and 4) from the fitting on both ends of the resin filter cartridges.
6. Replace both cartridges, taking care to insert them in the correct position and with the labels "FLOW" on the right-side (Fig. P).

NOTE

TO DISCONNECT THE PIPES FROM THE FITTINGS, PULL OUTWARDS WHILE KEEPING PRESSED THE FITTING LOCK RING.


Figure O


Figure P

7. Reconnect the pipes respecting the original positions, open the tap feeding the MILLDROP.
8. Modify the wire connection on the terminals board as described in the procedure for the purging of the hydraulic circuit (see Fig. F), direct the silicone pipe connected to the fitting C on the demineralizer back panel into the sink; insert the mains cable, turn on the equipment and let the water to drain for about 5 minutes.
9. Restore the electric connection on the electronic board and the lids.

Osmotic membrane

See Figure Q


NOTE

AFTER REPLACING THE OSMOTIC MEMBRANE AND THE PURGING OF THE HYDRAULIC CIRCUIT, REPLACE ALSO THE ACTIVE CARBON FILTER.

1. Turn off MILLDROP and remove the mains cable.
2. Proceed as with the resin filters and remove the whole filter assembly (A) from its bearings to access the osmotic membrane assembly (B).
3. Disconnect the pipes Ø8 (1) and (2) from the fittings and replace the osmotic membrane.
4. Reconnect the pipes respecting the original position.
5. Replace the active carbon filter as described on the next page.

NOTE

TO DISCONNECT THE PIPES FROM THE FITTINGS, PULL OUTWARDS WHILE KEEPING PRESSED THE FITTING LOCK RING.

**Figure Q**

6. Modify the wire connection on the terminals board as described in the procedure for the purging of the hydraulic circuit (see Fig. F), disconnect the pipe from the input resin filter fitting (see Fig. G) and perform the purging for 10 minutes.
7. Check for no leak on the equipment.
8. Restore the hydraulic and the electric connections, and replace the cover.
9. Check the quality of output demineralized water (see table).

Active carbon filter

See Figure R

NOTE

BESIDES THE ACTIVE CARBON FILTER, REPLACE ALSO THE OSMOTIC MEMBRANE.

1. Turn off MILLDROP and remove the mains cable.
2. Proceed as the resin filters replacement and extract the whole filters assemble from its bearings.
3. Disconnect the elbow Ø10 (1) and the "T" Ø10 (2) fitting from the carbon filter cartridge (A) and replace it.
4. Refer to Fig. Q on the previous page, and replace the osmotic membrane.
5. Reconnect the pipes respecting the original position.

NOTA

TO DISCONNECT THE PIPES FROM THE FITTINGS, PULL OUTWARDS WHILE KEEPING PRESSED THE FITTING LOCK RING.


Figure R

6. Modify the wire connection on the terminals board as described in the procedure for the purging of the hydraulic circuit (see Fig. F), disconnect the pipe from the input resin filter fitting (see Fig. G) and perform the purging for 10 minutes.
7. Restore the hydraulic and the electric connections, and replace the cover.
8. Check the quality of output demineralized water (see table).

REPLACING THE FLOW RESTRICTOR

See Figure S:

1. Turn off MILLDROP and remove the mains cable.
2. Turn off the tap feeding the water to MILLDROP, remove the lid and the board protection panel.
3. Lift and take out the filter assembly (A) from the brackets
4. Take out the osmotic membrane cartridge (B) and remove the red tube Ø6 from the “flow restrictor” union (C).
5. Unthread the flow restrictor and replace it with a new one (see Hydraulic diagram in Appendix D).

NOTE

REASSEMBLE THE FLOW RESTRICTOR UNION, AND SEAL USING “LOCTITE 5331” .

6. Restore the red tube Ø6 on the union, reassemble the osmotic membrane cartridge and the filter assembly.
7. Connect the mains cable and turn on MILLDROP.
8. Open the water feeding tap and check for no leaks.
9. Reassemble the board protection panel and the lid, fasten the screws.


Figure S

RECYCLING / DISPOSAL INSTRUCTIONS

Milldrop equipment is mainly built from fiber-reinforced polymers, metals and electronic components.

In case of disposal:

- Separate the various components according to the materials they are made of.
- Drop the equipment with a company specialized on recycling these products.
- Do not abandon the equipment in unsecured places.
- Always refer to current/applicable laws in the country of use.

INTRODUCTION

If you run into a problem or alarm while using the equipment, you should **not** be immediately concerned. It may not, in fact, be related to a breakdown but, more probably to an anomalous situation, often merely transitory (such as a blackout), or incorrect use.

In any case, it is important to first identify the cause of the anomaly and then take suitable corrective action, either autonomously or with the help of the **Technical Support Department** (see **Appendix Z**).

For this purpose, below, we provide instructions for diagnosing and resolving general problems.

At the starting on

It is possible that at the start of operation, and also in normal condition, the red LED (signaling poor outlet water quality) lights for a few seconds; this situation is normal.

On the contrary, if the LED remains lit for a time longer than 10 seconds (matched with the acoustic alarm), the equipment enters in blocking status; in this case it will be necessary to replace the resin filter and/or the osmotic membrane.

This fault condition is preceded by an intermediate phase during which Milldrop, even if correctly operating, turns on the yellow LED (matched with the acoustic alarm) to signal that the filter replacement is coming near.

Milldrop unused for a few days

Should the equipment be unused for a few days, the red LED could remain lit for a time longer than 10 seconds at the subsequent starting on, so resulting in the blocking of the equipment. In this case switch off the unit and switch on it again. If the problem occurs again, replace the filter.

NOTE

SHOULD YOU FORESEE NOT TO USE THE EQUIPMENT FOR MORE THAN 20 DAYS, WE SUGGEST TO EMPTY THE EXTERNAL TANK WHENEVER PROVIDED.

ANALYSIS AND RESOLUTION OF PROBLEMS


If your demineralizer is **not** working correctly, please make the following checks **before** calling the Technical Support Department:

PROBLEM	POSSIBLE CAUSE	PROPOSED SOLUTION
The demineralizer does not power-on.	The power cord is not plugged-in.	Plug it in.
	There is no voltage at the socket.	Check the cause of the lack of voltage at the socket and fix it.
	The main switch and/or differential switch are OFF.	Turn the switch ON.
	The mains fuses or the on-board fuse are blown.	Replace with new fuses of same rating value. (See <i>Summary Table</i> in Appendix A).
Alarm: water leak (red LED blinking)	Water leak inside the equipment.	Remove the water leak cause and dry the internal parts well.
	Electronic board failure.	Replace the electronic board.


PROBLEM	POSSIBLE CAUSE	PROPOSED SOLUTION
Alarm: lack of water (red LED blinking)	Entry tap closed.	Open the tap.
	Water work pressure insufficient.	Check the water pressure.
	Solenoid valve filter clogged up	Clean the filter (see rear view, page 5).
	Pressure switch failure.	Replace the pressure switch.
	Inlet solenoid valve failure.	Replace the solenoid valve.
	Active carbon filter clogged up.	Replace the active carbon filter.
Alarm: Poor outlet water quality (red LED blinking)	Resin filters depleted.	Replace the resin filter assembly.
	Osmotic membrane depleted.	Prepare to replace the osmotic membrane shortly
Alarm: Sufficient outlet water quality (yellow LED blinking)	Resin filters almost depleted.	Prepare to replace the resin filters shortly.
	Osmotic membrane almost depleted.	Prepare to replace the osmotic membrane shortly.
The demineralizer does not deliver any water.	Continuous working time exceeding 20 minutes.	Turn off the demineralizer and turn on again (or push again the tap's button).
	Start/stop signal cable detached	Connect the cable.
	Feeding tank's float failure.	Replace the float.
	Alarm signaling for inlet water lack or internal leak	Remove the problem.
	Lack of the power supply.	Remove the cause.
	Alarm signaling for poor outlet water quality.	Replace the resin filter assembly.
	Pump failure.	Replace the pump.
Low outlet flow rate	Pump failure.	Replace the pump.
	Resin filter and/or active carbon filter clogged up.	Replace the resin filter and/or the active carbon filter.
	Osmotic membrane clogged up.	Cleanse the osmotic membrane or replace it.

ELECTRIC DIAGRAM


Electronic board type 1

Electronic board type 2

HYDRAULIC DIAGRAM

Legend:

- | | | |
|--------------------------------------|---------------------------|------------------------|
| 1 – Water inlet solenoid valve | 4 – Pump assembly | 7 – Flow restrictor |
| 2 – Activated carbon filter assembly | 5 – Membrane assembly | 8 – Draining union |
| 3 – Pressure switch | 6 – Resin filter assembly | 9 – Conductivity probe |

DECLARATION OF CONFORMITY

Manufacturer's Name: **C.M.F. S.r.l. – Costruzioni Meccaniche Friulane**

Manufacturer's Address: **Via E.Ellero,9 – 33080 Cusano di Zoppola (PN [=Province of Pordenone]) – ITALY**

Distributor's Name: **M.O.COM. S.r.l. - Manifattura Odontoiatrica Complementare**

Distributor's Address: **Via delle Azalee, 1 - 20090 Buccinasco (MI) - ITALY**

Product description: **Reverse osmosis water demineralizer**

Model: 

Constructed in: **ITALY**

The undersigned declares that the above mentioned material


CONFORMS

to EEC Directives 89/336 – 73/23 (and subsequent updates).

Reference Standards : HD 195S5 EN 60555-2 EN 50081-1 EN 50082-1
EN 55022 EN 60204-1 CISPR 16-1 IEC 1000-4-2
IEC 1000-4-3 IEC 1000-4-4 EN 61010-1

02/01/2005

(Date)

Signature

Renzo Piccinin
Name and Surname

Legal Representative
Function

[illegible]

**FOR ANY REQUEST FOR
TECHNICAL SERVICE FOR THE PRODUCT,
WHETHER IN OR OUT OF WARRANTY,
DIRECTLY CONTACT THE

TECHNICAL SUPPORT DEPARTMENT

OF THE DEALER OR RESELLER
THAT SUPPLIED THE PRODUCT.**

M.O.COM. Srl is completely available to customers to provide any technical information about the product as well as to offer suggestions and advice on demineralization procedures.

In this regard, please refer to the following address:

M.O.COM. Srl
Customer Support
Via delle Azalee, 1
20090 Buccinasco (MI)
ITALY
Tel. (+39) 02-45701505
Fax (+39) 02-45701258
e-mail at@mocom.it
website www.mocom.it

To help us in the indispensable work of improving the quality of our products and service, please send your comments and/or suggestions by **letter** or **fax** or to the following **e-mail** address:

uc@mocom.it (Marketing/sale department)

Thank you in advance for your valuable assistance.

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